

# EGFR / HER1 / ErbB1 Antibody, Rabbit PAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 10001-RP01

## GENERAL INFORMATION

<b>Immunogen:</b>	Recombinant Human EGFR / ErbB / HER1 protein (Catalog#10001-H08H)
<b>Preparation</b>	Produced in rabbits immunized with purified, recombinant Human EGFR / ErbB / HER1 (rh EGFR / ErbB / HER1; Catalog#10001-H08H; Met 1-Ser 645; NP_005219). Total IgG was purified by Protein A affinity chromatography.
<b>Ig Type:</b>	Rabbit IgG
<b>Specificity:</b>	Human EGFR / ErbB / HER1
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## APPLICATIONS

<b>Applications:</b>	ELISA
----------------------	-------

## RECOMMENDED CONCENTRATION

<b>ELISA</b>	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human EGFR / HER1 / ErbB1.
--------------	---

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**