



Catalog Number: 10094-R004

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human GAPDH Protein (Catalog#10094-H07E)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human GAPDH (rh GAPDH; Catalog#10094-H07E; P04406; Met1-Glu335).
Ig Type:	Rabbit IgG
Clone ID:	004
Specificity:	Human GAPDH
Formulation:	Elution 3.0 +2.5% 2M Tris
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	
APPLICATIONS	
Applications:	WB
RECOMMENDED CONCENTRATION	
Western Blot	WB: 1:500-1:2000

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>



Catalog Number: 10094-R004

EliteRmab® is a registered trademark of Sino Biological Inc.

Anti-GAPDH rabbit monoclonal antibody at
1:500 dilution

Lane A: HeLa Whole Cell Lysate

Lane B: Jurkat Whole Cell Lysate

Lane C: HepG2 Whole Cell Lysate

Lane D: MCF7 Whole Cell Lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000
dilution.

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size:36 kDa

Observed band size:36 kDa

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>