

CSF1R / MCSF Receptor / CD115 Antibody (FITC), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10161-MM04-F

GENERAL INFORMATION

Immunogen:	Recombinant Human CSF1R / MCSF Receptor / CD115 Protein (Catalog#10161-H08H)
Reagents:	FITC-conjugated Mouse monoclonal antibody
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CSF1R / MCSF Receptor / CD115 (rh CSF1R / MCSF Receptor / CD115; Catalog#10161-H08H; NP_005202.2; Met1-Glu512) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.
Ig Type:	Mouse IgG1
Clone ID:	04
Specificity:	Human CSF1R / MCSF Receptor / CD115
Concentration:	10 µl/Test, 0.1 mg/ml
Formulation:	Aqueous solution containing 0.5% BSA and 0.09% sodium azide
Storage:	This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

APPLICATIONS

Applications: FCM

RECOMMENDED CONCENTRATION


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CSF1R / MCSF Receptor / CD115 Antibody (FITC), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10161-MM04-F



Flow cytometric analysis of Human M-CSFR(CD115) expression on human whole blood monocytes. Cells were stained with FITC-conjugated anti-Human M-CSFR(CD115). The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable monocytes.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.