

FAS / CD95 / APO-1 / TNFRSF6 Antibody, Rabbit MAb



Catalog Number: 10217-R225

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human CD95 / APO-1 / TNFRSF6 / FAS protein (Catalog#10217-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human CD95 / APO-1 / TNFRSF6 / FAS (rh CD95 / APO-1 / TNFRSF6 / FAS; Catalog#10217-H08H; NP_000034.1; Met1-Asn173).
Ig Type:	Rabbit IgG
Clone ID:	225
Specificity:	Human CD95 / APO-1 / TNFRSF6 / FAS
	No cross-reactivity in ELISA with Mouse FAS Rat FAS Human TNFR2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human FAS.

Please Note: Optimal concentrations/dilutions should be determined by the end user.