

**MMP9 / MMP-9 / CLG4B Antibody,
Rabbit MAb**

Catalog Number: 10327-R043

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human MMP9 / MMP-9 / CLG4B Protein (Catalog#10327-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human MMP9 / MMP-9 / CLG4B (rh MMP9 / MMP-9 / CLG4B; Catalog#10327-H08H; NP_004985.2; Met1-Asp707).
Ig Type:	Rabbit IgG
Clone ID:	043
Specificity:	Human MMP9 / MMP-9 / CLG4B
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA(Det)
RECOMMENDED CONCENTRATION	
Sandwich ELISA (Detection Ab)	ELISA(Det): 1:1000-1:10000 This antibody will detect Human MMP9/MMP-9/CLG4B in ELISA pair set (Catalog: # SEK10327 or Catalog: # SEKA10327). In a sandwich ELISA, it can be used as detection antibody when paired with (Catalog: # 10327-MM06).

Please Note: Optimal concentrations/dilutions should be determined by the end user.