

**MMP-9 / CLG4B Antibody, Rabbit MAb**

Catalog Number: 10327-RM04

EliteRmab® is a registered trademark of Sino Biological Inc.

| GENERAL INFORMATION       |   |
|---------------------------|---|
| <b>Immunogen:</b>         | Recombinant Human MMP9 protein (Catalog#10327-H08H)   |
| <b>Preparation</b>        | This antibody was obtained from a rabbit immunized with purified, recombinant Human matrix metalloproteinase 9 (rh MMP9; Catalog#10327-H08H; Met 1-Asp 707; NP_004985.2).   |
| <b>Ig Type:</b>           | Rabbit IgG  |
| <b>Clone ID:</b>          | 43  |
| <b>Specificity:</b>       | Human MMP-9 / CLG4B   |
| <b>Formulation:</b>       | 0.2 µm filtered solution in PBS   |
| <b>Storage:</b>           | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.<br>Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles. |
| APPLICATIONS              |   |
| <b>Applications:</b>      | ELISA   |
| RECOMMENDED CONCENTRATION |   |
| <b>ELISA</b>              | ELISA: 1:1000-1:2000<br>This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human MMP9.  |

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**