

MMP9/MMP-9/CLG4B Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10327-T08

GENERAL INFORMATION

Immunogen:	Recombinant Human MMP9/MMP-9/CLG4B Protein (Catalog#10327-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human MMP9/MMP-9/CLG4B (rh MMP9/MMP-9/CLG4B; Catalog#10327-H08H; NP_004985.2; Met1-Asp707). MMP9/MMP-9/CLG4B specific IgG was purified by Human MMP9/MMP-9/CLG4B affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human MMP9/MMP-9/CLG4B
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	IHC-P
----------------------	-------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:1000-1:5000
--------------	----------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

MMP9/MMP-9/CLG4B Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist

Catalog Number: 10327-T08



Immunochemical staining of human MMP9 in human tonsil with rabbit polyclonal antibody at 1:2000 dilution, formalin-fixed paraffin embedded sections.



Immunochemical staining of human MMP9 in human liver with rabbit polyclonal antibody at 1:2000 dilution, formalin-fixed paraffin embedded sections.