

CD50 / ICAM-3 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10333-MM07

GENERAL INFORMATION	
Immunogen:	10333-H08H
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human ICAM3 (rh ICAM3; Catalog#10333-H08H; CAA49369.1; Met 1-His 485). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	3B10B5B9
Specificity:	Human ICAM3 / CD50
	No cross-reactivity in ELISA with Human ICAM-2 / CD102 Human CD106 / VCAM1 Mouse CD106 / VCAM1 Mouse ICAM-1 / CD54
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	CD50
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000

Please Note: Optimal concentrations/dilutions should be determined by the end user.