

TRAIL R4 / CD264 / TNFRSF10D Antibody, Rabbit MAb



Catalog Number: 10413-R137

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human TRAIL R4 / CD264 / TNFRSF10D protein (Catalog#10413-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human TRAIL R4 / CD264 / TNFRSF10D (rh TRAIL R4 / CD264 / TNFRSF10D; Catalog#10413-H08H; NP_003831.2; Met1-His211).
Ig Type:	Rabbit IgG
Clone ID:	137
Specificity:	Human TRAIL R4 / CD264 / TNFRSF10D
	No cross-reactivity in ELISA with Human TNFRSF10A Human TNFRSF10B Human TNFSF10 Human TNF α Human TNF β
Formulation:	0.2 μ m filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human TRAIL R4 / CD264 / TNFRSF10D.

Please Note: Optimal concentrations/dilutions should be determined by the end user.