

# TRAIL R4 / CD264 / TNFRSF10D Antibody, Rabbit PAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 10413-RP01

GENERAL INFORMATION	
<b>Immunogen:</b>	Recombinant Human TRAIL R4 / CD264 / TNFRSF10D protein (Catalog#10413-H08H)
<b>Preparation</b>	Produced in rabbits immunized with purified, recombinant Human TRAIL R4 / CD264 / TNFRSF10D (rh TRAIL R4 / CD264 / TNFRSF10D; Catalog#10413-H08H; NP_003831.2; Met 1-His 211). Total IgG was purified by Protein A affinity chromatography.
<b>Ig Type:</b>	Rabbit IgG
<b>Specificity:</b>	Human TRAIL R4 / CD264 / TNFRSF10D
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
<b>Applications:</b>	ELISA
RECOMMENDED CONCENTRATION	
<b>ELISA</b>	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human TNFRSF10D.

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**