

# VE-Cadherin / CD144 / CDH5 Antibody, Rabbit PAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 10433-RP01

## GENERAL INFORMATION

<b>Immunogen:</b>	Recombinant Human VE-Cadherin / CD144 / CDH5 Protein
<b>Preparation</b>	Produced in rabbits immunized with purified, recombinant Human VE-Cadherin / CD144 / CDH5 (rh VE-Cadherin / CD144 / CDH5; CAA56306.1; 1M-593Q). Total IgG was purified by Protein A affinity chromatography.
<b>Ig Type:</b>	Rabbit IgG
<b>Specificity:</b>	Human VE-Cadherin / CD144 / CDH5
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## APPLICATIONS

**Applications:** ELISA

## RECOMMENDED CONCENTRATION

**ELISA** ELISA: 1:1000-1:2000  
This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human VE-Cadherin / CD144 / CDH5.

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**