

# Interferon Gamma/IFN gamma/IFNG Antibody, Rabbit PAb, Antigen Affinity Purified

Catalog Number: 105995-T08



Sino Biological  
Biological Solution Specialist

## GENERAL INFORMATION

<b>Immunogen:</b>	A synthetic peptide corresponding to the C-terminus of the Mouse Interferon Gamma/IFN gamma/IFNG
<b>Preparation</b>	Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Mouse Interferon Gamma/IFN gamma/IFNG, and purified by antigen affinity chromatography.
<b>Ig Type:</b>	Rabbit IgG
<b>Specificity:</b>	Mouse Interferon Gamma/IFN gamma/IFNG
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## APPLICATIONS

<b>Applications:</b>	IHC-P
----------------------	-------

## RECOMMENDED CONCENTRATION

<b>IHC-P</b>	IHC-P: 1:200-1:1000
--------------	---------------------


***Please Note: Optimal concentrations/dilutions should be determined by the end user.***

**Interferon Gamma/IFN gamma/IFNG  
Antibody, Rabbit PAb, Antigen  
Affinity Purified**


Catalog Number: 105995-T08



Sino Biological  
Biological Solution Specialist



Immunohistochemical staining of mouse IFNG in mouse kidney with rabbit polyclonal antibody at 1:500 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of mouse IFNG in mouse testis with rabbit polyclonal antibody at 1:500 dilution, formalin-fixed paraffin embedded sections.