

# TNF-alpha / TNFA / TNFSF2 Antibody (Biotin), Rabbit MAb



Catalog Number: 10602-RE23-B

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	10602-HNAE
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human TNF-alpha / TNFA (rh TNF-alpha / TNFA; Catalog#10602-HNAE; NP_000585.2; Val77-Leu233) and then biotinylated.
Ig Type:	Rabbit IgG
Clone ID:	E23
Specificity:	Human TNF-alpha / TNFA
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	TNF
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 0.1-1 µg/ml

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**