

CSF2RA / GM-CSFR / CD116 Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10701-RP01

GENERAL INFORMATION

Immunogen:	10701-H08H
Preparation	Produced in rabbits immunized with purified, recombinant Human CD116 / GM-CSFR (rh GM-CSFR; Catalog#10701-H08H; P15509-1; Met 1-Gly 320). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human CD116 / GM-CSFR
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	CD116

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000
--------------	----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.