

**N-Cadherin / CD325 / CDH2 Antibody (FITC), Rabbit MAb**

Catalog Number: 11039-R020-F

EliteRmab® is a registered trademark of Sino Biological Inc.

**GENERAL INFORMATION**

<b>Immunogen:</b>	Recombinant Human N-Cadherin / CD325 / CDH2 protein (Catalog#11039-H08H)
<b>Reagents:</b>	FITC-conjugated Rabbit monoclonal antibody
<b>Preparation</b>	This antibody was obtained from a rabbit immunized with purified, recombinant Human N-Cadherin / CD325 / CDH2 (rh N-Cadherin / CD325 / CDH2; Catalog#11039-H08H; NP_001783.2; Met1-Ala724) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.
<b>Ig Type:</b>	Rabbit IgG
<b>Clone ID:</b>	020
<b>Specificity:</b>	Human N-Cadherin / CD325 / CDH2
<b>Concentration:</b>	10 µl/Test, 0.1 mg/ml
<b>Formulation:</b>	PBS solution containing 0.5% BSA and 0.09% sodium azide
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
<b>Alternative Names:</b>	CDH2

**APPLICATIONS**

<b>Applications:</b>	FCM
----------------------	-----


**RECOMMENDED CONCENTRATION***Please Note: Optimal concentrations/dilutions should be determined by the end user.*

# N-Cadherin / CD325 / CDH2 Antibody (FITC), Rabbit MAb



Catalog Number: 11039-R020-F

EliteRmab® is a registered trademark of Sino Biological Inc.



Flow cytometry analysis of N-Cadherin on HeLa cells. HeLa cells were harvested with (Left) or without (Right) trypsinization (please note, the epitope is sensitive to trypsin) and stained with either Rabbit IgG isotype control FITC (dashed line) or Rabbit Anti-Human CD325 antibody FITC (solid line) at matching concentrations.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to [www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html](http://www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html) for technical protocols.