

CD16/Fc gamma RIII Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11046-MM10

GENERAL INFORMATION

Immunogen:	Recombinant Human CD16b/Fc gamma RIIIb Protein (Catalog#11046-H08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD16b/Fc gamma RIIIb (rh CD16b/Fc gamma RIIIb; Catalog#11046-H08H; O75015.2; Met1-Ser200). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	10
Specificity:	Human CD16 (CD16a & CD16b)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION

Flow Cytometry	FCM: 1:25-1:100
-----------------------	-----------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CD16/Fc gamma RIII Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11046-MM10



Flow cytometric analysis of Human CD16 expression on human whole blood lymphocytes. Cells were stained with purified anti-Human CD16, then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.