

# Influenza A H5N1 (A/Anhui/1/2005) Hemagglutinin / HA Antibody, Rabbit PAb, Antigen Affinity Purified

Catalog Number: 11048-T30



Sino Biological  
Biological Solution Specialist

## GENERAL INFORMATION

<b>Immunogen:</b>	Recombinant Influenza A H5N1 (A/Anhui/1/2005) Hemagglutinin / HA Protein (Catalog#11048-V08H1)
<b>Preparation</b>	Produced in rabbits immunized with purified, recombinant Influenza A H5N1 (A/Anhui/1/2005) Hemagglutinin / HA ( Catalog#11048-V08H1; ABD28180.1; Met1-Gln527). Influenza A H5N1 (A/Anhui/1/2005) Hemagglutinin / HA specific IgG was purified by Influenza A H5N1 (A/Anhui/1/2005) Hemagglutinin / HA affinity chromatography.
<b>Ig Type:</b>	Rabbit IgG
<b>Specificity:</b>	Influenza A H5N1 (A/Anhui/1/2005) Hemagglutinin / HA
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

## APPLICATIONS

<b>Applications:</b>	ELISA, IHC-P, FCM, ICC/IF, IP
	(Antibody's applications have not been validated with corresponding viruses. Optimal concentrations/dilutions should be determined by the end user.)

## RECOMMENDED CONCENTRATION

<b>ELISA</b>	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Influenza A H5N1 (A/Anhui/1/2005) Hemagglutinin / HA.
--------------	--

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**

**For Research Use Only. Not for use in diagnostic or therapeutic procedures.**

**Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)**

**Website: <http://www.sinobiological.com>**