

Influenza H5N1 (A/turkey/Turkey/1/2005) Hemagglutinin / HA Antibody, Rabbit

Product Number: 11061-RP01



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION	
Immunogen:	Recombinant Influenza H5N1 (A/turkey/Turkey/1/2005) Hemagglutinin / HA Protein (Catalog#11061-V08H1)
Preparation	Produced in rabbits immunized with purified, recombinant Influenza H5N1 (A/turkey/Turkey/1/2005) Hemagglutinin / HA (Catalog#11061-V08H1; ABD73284.1; Met1-Gln531). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Influenza H5N1 (A/turkey/Turkey/1/2005) Hemagglutinin / HA
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	Hemagglutinin,HA
APPLICATIONS	
Applications:	ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Influenza H5N1 (A/turkey/Turkey/1/2005) Hemagglutinin / HA.

Please Note: Optimal concentrations/dilutions should be determined by the end user.