

# Influenza A H1N1 (A/Brevig Mission/1/1918) Hemagglutinin / HA Antibody, Rabbit PAb

Catalog Number: 11068-RP01



Sino Biological  
Biological Solution Specialist

## GENERAL INFORMATION

<b>Immunogen:</b>	Influenza A H1N1 (A/Brevig Mission/1/1918) Hemagglutinin / HA protein (Catalog#11068-V08H)
<b>Preparation</b>	Produced in rabbits immunized with purified, recombinant Influenza A H1N1 (A/Brevig Mission/1/1918) Hemagglutinin / HA (Catalog#11068-V08H; AAD17229.1; Met1-Gln529). Total IgG was purified by Protein A affinity chromatography.
<b>Ig Type:</b>	Rabbit IgG
<b>Specificity:</b>	Influenza A H1N1 (A/Brevig Mission/1/1918) Hemagglutinin / HA
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
<b>Alternative Names:</b>	Hemagglutinin,HA

## APPLICATIONS

<b>Applications:</b>	ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)

## RECOMMENDED CONCENTRATION

<b>ELISA</b>	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Influenza A H1N1 (A/Brevig Mission/1/1918) Hemagglutinin / HA.
--------------	---

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**