

Anti-H7N9 Hemagglutinin / HA Antibody (Biotin), Hemagglutination Inhibition Activity, Rabbit MAb

Catalog Number: 11082-R002-B



Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Ig Type:	Rabbit IgG
Clone ID:	002
Specificity:	Anti-H7N9 HA / Hemagglutinin
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	HA
APPLICATIONS	
Applications:	ELISA
	(Antibody's applications have not been validated with corresponding viruses. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 0.1-1 µg/ml

Please Note: Optimal concentrations/dilutions should be determined by the end user.