

IGSF3 Antibody (APC), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11290-MM05-A

GENERAL INFORMATION

Immunogen:	Recombinant Human IGSF3 protein (Catalog#11290-H08H)
Reagents:	APC-conjugated Mouse monoclonal antibody
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human IGSF3 (rh IGSF3; Catalog#11290-H08H; O75054-1; Met1-Ala1124) and conjugated with APC under optimum conditions, the unreacted APC was removed.
Ig Type:	Mouse IgG2a
Clone ID:	05
Specificity:	Human IGSF3
Concentration:	10 µl/Test, 0.1 mg/ml
Formulation:	Aqueous solution containing 0.5% BSA and 0.09% sodium azide
Storage:	This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION


Please Note: Optimal concentrations/dilutions should be determined by the end user.

IGSF3 Antibody (APC), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11290-MM05-A



Profile of anti-Human IGSF3 reactivity on A549 cells analyzed by flow cytometry.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system.

Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.