

IGSF3 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 11290-T08

GENERAL INFORMATION

Immunogen:	Recombinant Human IGSF3 Protein (Catalog#11290-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human IGSF3 (rh IGSF3; Catalog#11290-H08H; O75054-1; Met1-Ala1124). IGSF3 specific IgG was purified by Human IGSF3 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human IGSF3
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications: IHC-P

RECOMMENDED CONCENTRATION

IHC-P IHC-P: 1:500-1:2000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

IGSF3 Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist

Catalog Number: 11290-T08



Immunochemical staining of human IGSF3 in human kidney with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections).



IGSF3 Antibody

IGSF3 Antigen

Immunochemical staining of human IGSF3 in human brain with rabbit polyclonal antibody (1:1000, formalin-fixed paraffin embedded sections). The left panel: tissue incubated with primary antibody; The right panel: tissue incubated with the mixture of primary