

Influenza A H1N1 (A/New Caledonia/20/1999) Hemagglutinin / HA Antibody, Rabbit PAb

Catalog Number: 11683-RP01



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION

Immunogen:	Recombinant Influenza A H1N1 (A/New Caledonia/20/1999) Hemagglutinin / HA Protein (Catalog#11683-V08H)
Preparation	Produced in rabbits immunized with purified, recombinant Influenza A H1N1 (A/New Caledonia/20/1999) Hemagglutinin / HA (Catalog#11683-V08H; AAP34324.1; Met1-Gln528). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Influenza A H1N1 (A/New Caledonia/20/1999) Hemagglutinin / HA
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	Hemagglutinin,HA

APPLICATIONS

Applications:	ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Influenza A H1N1 (A/New Caledonia/20/1999) Hemagglutinin / HA.
--------------	---

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>