

Influenza H1N1 (A/Puerto Rico/8/1934) Hemagglutinin / HA Antibody, Rabbit MAb

Catalog Number: 11684-R016



Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Influenza H1N1 (A/Puerto Rico/8/1934) Hemagglutinin / HA Protein (Catalog#11684-V08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Influenza H1N1 (A/Puerto Rico/8/1934) Hemagglutinin / HA (Catalog#11684-V08H; ABD77675.1; Met1-Gln528).
Ig Type:	Rabbit IgG
Clone ID:	016
Specificity:	Influenza H1N1 (A/Puerto Rico/8/1934) Hemagglutinin / HA
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA(Cap)
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
Sandwich ELISA (Capture Ab)	ELISA(Cap): 1:250-1:2000 This antibody will detect Influenza H1N1 (A/Puerto Rico/8/1934) Hemagglutinin / HA in ELISA pair set (Catalog: # SEKA11684). In a sandwich ELISA, it can be used as capture antibody when paired with (Catalog: # 11684-R107).

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>