

Influenza A H1N1 (A/Puerto Rico/8/1934) Hemagglutinin / HA Antibody, Rabbit PAb, Antigen Affinity Purified

Catalog Number: 11684-T62



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION

Immunogen:	Recombinant Influenza A H1N1 (A/Puerto Rico/8/1934) Hemagglutinin / HA Protein (Catalog#11684-V08H)
Preparation	Produced in rabbits immunized with purified, recombinant Influenza A H1N1 (A/Puerto Rico/8/1934) Hemagglutinin / HA (Catalog#11684-V08H; ABD77675.1; Met1-Gln528). Influenza A H1N1 (A/Puerto Rico/8/1934) Hemagglutinin / HA specific IgG was purified by Influenza A H1N1 (A/Puerto Rico/8/1934) Hemagglutinin / HA affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Influenza A H1N1 (A/Puerto Rico/8/1934) Hemagglutinin / HA
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	Hemagglutinin,HA

APPLICATIONS

Applications:	WB,ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:2000-1:5000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Influenza A H1N1 (A/Puerto Rico/8/1934) Hemagglutinin / HA.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

**Influenza A H1N1 (A/Puerto Rico/8/1934)
Hemagglutinin / HA Antibody, Rabbit
PAb, Antigen Affinity Purified**

Catalog Number: 11684-T62



Sino Biological
Biological Solution Specialist



Anti-Influenza A H1N1 (A/Puerto Rico/8/1934)
Hemagglutinin / HA rabbit polyclonal antibody
at 1:2000 dilution
Sample: Influenza A H1N1 (A/Puerto
Rico/8/1934) Hemagglutinin / HA Recombinant
Protein
Lane A: 20ng
Lane B: 5ng

Secondary
Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000
dilution.

Developed using the ECL technique.
Performed under reducing conditions.