

**Influenza H5N1 (A/Common magpie/Hong Kong/2256/2006) Hemagglutinin / HA Antibody, Rabbit PAb**



**Sino Biological**  
Biological Solution Specialist

Catalog Number: 11700-RP01

**GENERAL INFORMATION**

<b>Immunogen:</b>	Recombinant Influenza H5N1 (A/Common magpie/Hong Kong/2256/2006) Hemagglutinin / HA Protein (Catalog#11700-V08H)
<b>Preparation</b>	Produced in rabbits immunized with purified, recombinant Influenza H5N1 (A/Common magpie/Hong Kong/2256/2006) Hemagglutinin / HA ( Catalog#11700-V08H; ABJ96777.1; Met1-Glu530). Total IgG was purified by Protein A affinity chromatography.
<b>Ig Type:</b>	Rabbit IgG
<b>Specificity:</b>	Influenza H5N1 (A/Common magpie/Hong Kong/2256/2006) Hemagglutinin / HA
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
<b>Alternative Names:</b>	Hemagglutinin,HA

**APPLICATIONS**

<b>Applications:</b>	ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)

**RECOMMENDED CONCENTRATION**

<b>ELISA</b>	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Influenza H5N1 (A/Common magpie/Hong Kong/2256/2006) Hemagglutinin / HA.
--------------	---

***Please Note: Optimal concentrations/dilutions should be determined by the end user.***