

# IFN gamma / IFNG Antibody, Rabbit MAb



Catalog Number: 11725-R003

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
<b>Immunogen:</b>	Recombinant Human IFN gamma / IFNG Protein (Catalog#11725-HNAS)
<b>Preparation</b>	This antibody was obtained from a rabbit immunized with purified, recombinant Human IFN gamma / IFNG (rh IFN gamma / IFNG; Catalog#11725-HNAS; NP_000610.2; Met1-Gln166).
<b>Ig Type:</b>	Rabbit IgG
<b>Clone ID:</b>	003
<b>Specificity:</b>	Human IFN gamma / IFNG
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
<b>Applications:</b>	FCM
RECOMMENDED CONCENTRATION	
<b>Flow Cytometry</b>	FCM: 1:25-1:100

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**


# IFN gamma / IFNG Antibody, Rabbit MAb

Catalog Number: 11725-R003



Sino Biological  
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.



Flow cytometric analysis of Human IFN- $\gamma$  expression on human peripheral blood lymphocytes. Human peripheral blood mononuclear cells were stimulated for 4-6 hours with PMA and Ionomycin in the presence of GolgiPlug. The cells were treated according to manufacturer's manual (BD Pharmingen™ Cat. No. 554714), stained with purified anti-Human IFN- $\gamma$ , then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to [www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html](http://www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html) for technical protocols.