

# Interferon Gamma/IFN gamma/IFNG Antibody (Biotin), Rabbit MAb



Catalog Number: 11725-R209-B

EliteRmab® is a registered trademark of Sino Biological Inc.

| GENERAL INFORMATION       |  |
|---------------------------|--|
| <b>Immunogen:</b>         | Recombinant Human Interferon Gamma/IFN gamma/IFNG Protein (Catalog#11725-HNAS)   |
| <b>Preparation</b>        | This antibody was obtained from a rabbit immunized with purified, recombinant Human Interferon Gamma/IFN gamma/IFNG (rh Interferon Gamma/IFN gamma/IFNG; Catalog#11725-HNAS; NP_000610.2; Met1-Gln166) and then biotinylated.                      |
| <b>Ig Type:</b>           | Rabbit IgG   |
| <b>Clone ID:</b>          | 209  |
| <b>Specificity:</b>       | Human Interferon Gamma/IFN gamma/IFNG  |
| <b>Formulation:</b>       | 0.2 µm filtered solution in PBS  |
| <b>Storage:</b>           | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles. |
| <b>Alternative Names:</b> | IFG, IFI, IFN gamma, Interferon Gamma  |
| APPLICATIONS              |  |
| <b>Applications:</b>      | ELISA  |
| RECOMMENDED CONCENTRATION |  |
| <b>ELISA</b>              | ELISA: 0.1-1 µg/ml   |

*Please Note: Optimal concentrations/dilutions should be determined by the end user.*