

# Interferon Gamma/IFN gamma/IFNG Antibody (HRP), Rabbit MAb



Catalog Number: 11725-R209-H

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
<b>Immunogen:</b>	Recombinant Human Interferon Gamma/IFN gamma/IFNG Protein (Catalog#11725-HNAS)
<b>Preparation</b>	This antibody was obtained from a rabbit immunized with purified, recombinant Human Interferon Gamma/IFN gamma/IFNG (rh Interferon Gamma/IFN gamma/IFNG; Catalog#11725-HNAS; NP_000610.2; Met1-Gln166), and conjugated with horseradish-peroxidase (HRP).
<b>Ig Type:</b>	Rabbit IgG
<b>Clone ID:</b>	209
<b>Specificity:</b>	Human Interferon Gamma/IFN gamma/IFNG
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
<b>Applications:</b>	ELISA
RECOMMENDED CONCENTRATION	
<b>ELISA</b>	ELISA: 0.1-1 µg/ml

*Please Note: Optimal concentrations/dilutions should be determined by the end user.*