

# Interferon Gamma/IFN gamma/IFNG Antibody, Rabbit MAb



Catalog Number: 11725-R209

EliteRmab® is a registered trademark of Sino Biological Inc.

## GENERAL INFORMATION

|                     |  |
|---------------------|--|
| <b>Immunogen:</b>   | Recombinant Human Interferon Gamma/IFN gamma/IFNG Protein (Catalog#11725-HNAS)   |
| <b>Preparation</b>  | This antibody was obtained from a rabbit immunized with purified, recombinant Human Interferon Gamma/IFN gamma/IFNG (rh Interferon Gamma/IFN gamma/IFNG; Catalog#11725-HNAS; NP_000610.2; Met1-Gln166).  |
| <b>Ig Type:</b>     | Rabbit IgG   |
| <b>Clone ID:</b>    | 209  |
| <b>Specificity:</b> | Human Interferon Gamma/IFN gamma/IFNG  |
| <b>Formulation:</b> | 0.2 µm filtered solution in PBS  |
| <b>Storage:</b>     | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles. |

## APPLICATIONS

|                      |            |
|----------------------|------------|
| <b>Applications:</b> | ELISA(Det) |
|----------------------|------------|

## RECOMMENDED CONCENTRATION

|                                      |   |
|--------------------------------------|---|
| <b>Sandwich ELISA (Detection Ab)</b> | ELISA(Det): 1:1000-1:10000<br>This antibody will detect Human Interferon Gamma/IFN gamma/IFNG in ELISA pair set (Catalog: # SEK11725A). In a sandwich ELISA, it can be used as detection antibody when paired with (Catalog: # 11725-R238). |
|--------------------------------------|---|

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**