

# Interferon Gamma/IFN gamma/IFNG Antibody, Rabbit MAb



Catalog Number: 11725-R238

EliteRmab® is a registered trademark of Sino Biological Inc.

## GENERAL INFORMATION

<b>Immunogen:</b>	Recombinant Human Interferon Gamma/IFN gamma/IFNG Protein (Catalog#11725-HNAS)
<b>Preparation</b>	This antibody was obtained from a rabbit immunized with purified, recombinant Human Interferon Gamma/IFN gamma/IFNG (rh Interferon Gamma/IFN gamma/IFNG; Catalog#11725-HNAS; NP_000610.2; Met1-Gln166).
<b>Ig Type:</b>	Rabbit IgG
<b>Clone ID:</b>	238
<b>Specificity:</b>	Human Interferon Gamma/IFN gamma/IFNG
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

## APPLICATIONS

<b>Applications:</b>	ELISA(Cap)
----------------------	------------

## RECOMMENDED CONCENTRATION

<b>Sandwich ELISA (Capture Ab)</b>	ELISA(Cap): 1:250-1:2000 This antibody will detect Human Interferon Gamma/IFN gamma/IFNG in ELISA pair set (Catalog: # SEK11725A). In a sandwich ELISA, it can be used as capture antibody when paired with (Catalog: # 11725-R209).
--	---

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**