

CD62L / L-Selectin / SELL Antibody (HRP), Rabbit MAb



Catalog Number: 11838-R018-H

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human CD62L / L-Selectin / SELL protein (Catalog#11838-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human CD62L / L-Selectin / SELL (rh CD62L / L-Selectin / SELL; Catalog#11838-H08H; P14151-1; Met 1-Asn 332), and conjugated with horseradish-peroxidase (HRP).
Ig Type:	Rabbit IgG
Clone ID:	018
Specificity:	Human CD62L / L-Selectin / SELL
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	SELL
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 0.1-1 µg/ml

Please Note: Optimal concentrations/dilutions should be determined by the end user.