

CD62L / L-Selectin Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11838-RP01

GENERAL INFORMATION

| | |
|---------------------|---|
| Immunogen: | Recombinant Human CD62L / L-Selectin protein (Catalog#11838-H08H) |
| Preparation | Produced in rabbits immunized with purified, recombinant Human CD62L / L-Selectin (rh SELL; Catalog#11838-H08H; P14151-1; Met1-Asn332). Total IgG was purified by Protein A affinity chromatography. |
| Ig Type: | Rabbit IgG |
| Specificity: | Human CD62L / L-Selectin |
| Formulation: | 0.2 µm filtered solution in PBS |
| Storage: | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles. |

APPLICATIONS

| | |
|----------------------|-------|
| Applications: | ELISA |
|----------------------|-------|

RECOMMENDED CONCENTRATION

| | |
|--------------|--|
| ELISA | ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human CD62L / L-Selectin. |
|--------------|--|

Please Note: Optimal concentrations/dilutions should be determined by the end user.