

# IL-2 Antibody (HRP), Rabbit PAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 11848-RP01-H

## GENERAL INFORMATION

<b>Immunogen:</b>	Recombinant Human IL-2 protein (Catalog#11848-HNAE)
<b>Preparation</b>	Produced in rabbits immunized with purified, recombinant Human IL-2 (rh IL-2; Catalog#11848-HNAE; P60568-1; Ala21-Thr153). Total IgG was purified by Protein A affinity chromatography and conjugated with horseradish-peroxidase (HRP).
<b>Ig Type:</b>	Rabbit IgG
<b>Specificity:</b>	Human IL-2
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
<b>Alternative Names:</b>	IL-2, Interleukin-2, lymphokine, TCGF

## APPLICATIONS

**Applications:** ELISA

## RECOMMENDED CONCENTRATION

**ELISA** ELISA: 0.1-1 µg/ml

*Please Note: Optimal concentrations/dilutions should be determined by the end user.*