

IL2/IL-2/Interleukin-2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 11848-T16

GENERAL INFORMATION

Immunogen:	Recombinant Human IL2/IL-2/Interleukin-2 Protein (Catalog#11848-HNAY1)
Preparation	Produced in rabbits immunized with purified, recombinant Human IL2/IL-2/Interleukin-2 (rh IL2/IL-2/Interleukin-2; Catalog#11848-HNAY1; NP_000577.2; Ala 21-Thr153, Leu 38 Met, Leu 39 Ser, Cys 145 Ala). IL2/IL-2/Interleukin-2 specific IgG was purified by Human IL2/IL-2/Interleukin-2 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human IL2/IL-2/Interleukin-2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications: ELISA

RECOMMENDED CONCENTRATION

ELISA ELISA: 1:5000-1:10000
This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human IL2/IL-2/Interleukin-2.

Please Note: Optimal concentrations/dilutions should be determined by the end user.