

Human coronavirus spike glycoprotein Antibody, Rabbit MAb



Catalog Number: 40021-R028

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human coronavirus Human coronavirus spike glycoprotein protein (Catalog#40021-V08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human coronavirus Human coronavirus spike glycoprotein (Catalog#40021-V08H; YP_173238.1; Met1-Arg760).
Ig Type:	Rabbit IgG
Clone ID:	028
Specificity:	Human coronavirus Human coronavirus spike glycoprotein
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	S
APPLICATIONS	
Applications:	ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human coronavirus Human coronavirus spike glycoprotein.

Please Note: Optimal concentrations/dilutions should be determined by the end user.