

MERS-CoV Nucleoprotein / NP Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 40068-RP01

GENERAL INFORMATION

Immunogen:	Recombinant MERS-CoV (NCoV / Novel coronavirus) Nucleoprotein / NP protein (Catalog#40068-V08B)
Preparation	Produced in rabbits immunized with purified, recombinant MERS-CoV (NCoV / Novel coronavirus) Nucleoprotein / NP (Catalog#40068-V08B; AFS88943.1; Met1-Asp413). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	MERS-CoV (NCoV / Novel coronavirus) Nucleoprotein / NP
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	NP

APPLICATIONS

Applications:	ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect MERS-CoV (NCoV / Novel coronavirus) Nucleoprotein / NP.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.