

Influenza A H7N9 (A/Anhui/1/2013) Hemagglutinin / HA Antibody, Mouse MAb

Catalog Number: 40103-MM06



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION

Immunogen:	Recombinant Influenza A H7N9 (A/Anhui/1/2013) Hemagglutinin / HA Protein (Catalog#40103-V08B)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant H7N9 HA (Catalog#40103-V08B; EP1439507; Met1-Val524). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	06
Specificity:	Influenza A H7N9 (A/Shanghai/1/2013) Hemagglutinin / HA No cross-reactivity in ELISA with Human cell lysate (293 cell line)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	Hemagglutinin,HA

APPLICATIONS

Applications:	ELISA IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)
----------------------	---

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Influenza A H7N9 (A/Shanghai/1/2013) Hemagglutinin / HA.
--------------	---

Please Note: Optimal concentrations/dilutions should be determined by the end user.