

# SARS-CoV / SARS-CoV-2 Nucleocapsid Antibody (HRP), Rabbit MAb

Catalog Number: 40143-R004-H



Sino Biological  
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
<b>Immunogen:</b>	40143-V08B
<b>Preparation</b>	This product is a recombinant monoclonal antibody expressed from HEK293 cells and conjugated with horseradish-peroxidase (HRP).
<b>Ig Type:</b>	Rabbit IgG
<b>Clone ID:</b>	004
<b>Specificity:</b>	SARS-CoV Nucleocapsid Protein
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
<b>Alternative Names:</b>	NP-CoV
APPLICATIONS	
<b>Applications:</b>	ELISA
	(Antibody's applications have not been validated with corresponding viruses. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
<b>ELISA</b>	ELISA: 0.1-1 µg/ml

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**