

# SARS-CoV / SARS-CoV-2 Spike Antibody (Biotin), Rabbit MAb



Catalog Number: 40150-R007-B

EliteRmab® is a registered trademark of Sino Biological Inc.

| GENERAL INFORMATION       |  |
|---------------------------|--|
| <b>Immunogen:</b>         | 40150-V08B1  |
| <b>Preparation</b>        | This product is a recombinant monoclonal antibody expressed from HEK293 cells and then biotinylated.   |
| <b>Ig Type:</b>           | Rabbit IgG   |
| <b>Clone ID:</b>          | 007  |
| <b>Specificity:</b>       | SARS-CoV-2 (2019-nCoV) Spike S1 Protein<br>SARS-CoV-2 (2019-nCoV) Spike RBD Protein<br>SARS-CoV Spike S1 Protein<br>SARS-CoV Spike RBD Protein   |
| <b>Formulation:</b>       | 0.2 µm filtered solution in PBS  |
| <b>Storage:</b>           | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles. |
| <b>Alternative Names:</b> | Spike  |
| APPLICATIONS              |  |
| <b>Applications:</b>      | ELISA<br><br>(Antibody's applications have not been validated with corresponding viruses. Optimal concentrations/dilutions should be determined by the end user.)  |
| RECOMMENDED CONCENTRATION |  |
| <b>ELISA</b>              | ELISA: 0.1-1 µg/ml   |

*Please Note: Optimal concentrations/dilutions should be determined by the end user.*