

SARS-CoV / SARS-CoV-2 Spike Antibody (HRP), Rabbit MAb



Catalog Number: 40150-R007-H

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	40150-V08B1
Preparation	This product is a recombinant monoclonal antibody expressed from HEK293 cells and conjugated with horseradish-peroxidase (HRP).
Ig Type:	Rabbit IgG
Clone ID:	007
Specificity:	SARS-CoV-2 (2019-nCoV) Spike S1 Protein SARS-CoV-2 (2019-nCoV) Spike RBD Protein SARS-CoV Spike S1 Protein SARS-CoV Spike RBD Protein
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	Spike
APPLICATIONS	
Applications:	ELISA
	(Antibody's applications have not been validated with corresponding viruses. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 0.1-1 µg/ml

Please Note: Optimal concentrations/dilutions should be determined by the end user.