

Spike Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 40589-T30

GENERAL INFORMATION

Immunogen:	Recombinant SARS-CoV-2 (2019-nCoV) Spike Protein (Catalog#40589-V08B1)
Preparation	Produced in rabbits immunized with purified, recombinant SARS-CoV-2 (2019-nCoV) Spike (Catalog#40589-V08B1; YP_009724390.1; Val 16-Pro1213). Spike specific IgG was purified by SARS-CoV-2 (2019-nCoV) Spike affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	SARS-CoV-2 (2019-nCoV) Spike
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect SARS-CoV-2 (2019-nCoV) Spike.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>