

SARS-CoV-2 Spike Antibody (HRP), Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 40591-T62-H

GENERAL INFORMATION

Immunogen:	40591-V08H
Preparation	Produced in rabbits immunized with purified, recombinant SARS-CoV-2 (2019-nCoV) Spike S1-His Recombinant Protein (HPLC-verified) (Catalog#40591-V08H; YP_009724390.1; Val16-Arg685). The specific IgG was purified by SARS-CoV-2 / 2019-nCoV Spike affinity chromatography and conjugated with horseradish-peroxidase (HRP).
Ig Type:	Rabbit IgG
Specificity:	SARS-CoV-2 (2019-nCoV) Spike S-His Protein (Catalog#40591-V08H)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	Spike

APPLICATIONS

Applications:	ELISA
	(Antibody's applications have not been validated with corresponding viruses. Optimal concentrations/dilutions should be determined by the end user.)

RECOMMENDED CONCENTRATION

ELISA	ELISA: 0.1-1 µg/ml
--------------	--------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.