

SARS-CoV-2 (2019-nCoV) Spike RBD Antibody, Mouse PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 40592-MP01

GENERAL INFORMATION

Immunogen:	Recombinant SARS-CoV-2 (2019-nCoV) Spike RBD-His Protein (Catalog#40592-V08H)
Preparation	Produced in mice immunized with purified, recombinant SARS-CoV-2 / 2019-nCoV Spike/RBD Protein (Catalog#40592-V08H; YP_009724390.1; Arg319-Phe541).
Specificity:	SARS-CoV-2 (2019-nCoV) Spike RBD (His Tag) (Cat# 40592-V08H)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	Spike

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect SARS-CoV-2 (2019-nCoV) Spike RBD Protein.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.