

# FAS / CD95 / APO-1 / TNFRSF6 Antibody (APC), Rabbit MAb



Catalog Number: 50027-R004-A

EliteRmab® is a registered trademark of Sino Biological Inc.

## GENERAL INFORMATION

<b>Immunogen:</b>	Recombinant Mouse CD95 / APO-1 / TNFRSF6 / FAS protein (Catalog#50027-M08H)
<b>Reagents:</b>	APC-conjugated Rabbit monoclonal antibody
<b>Preparation</b>	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD95 / APO-1 / TNFRSF6 / FAS (rM CD95 / APO-1 / TNFRSF6 / FAS; Catalog#50027-M08H; NP_032013.2; Met1-Arg169) and conjugated with APC under optimum conditions, the unreacted APC was removed.
<b>Ig Type:</b>	Rabbit IgG
<b>Clone ID:</b>	004
<b>Specificity:</b>	Mouse CD95 / APO-1 / TNFRSF6 / FAS
<b>Concentration:</b>	5 µl/Test, 0.1 mg/ml
<b>Formulation:</b>	PBS solution containing 0.5% BSA and 0.09% sodium azide
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
<b>Alternative Names:</b>	FAS

## APPLICATIONS

<b>Applications:</b>	FCM
----------------------	-----


## RECOMMENDED CONCENTRATION

*Please Note: Optimal concentrations/dilutions should be determined by the end user.*

**FAS / CD95 / APO-1 / TNFRSF6  
Antibody (APC), Rabbit MAb**

Catalog Number: 50027-R004-A

EliteRmab® is a registered trademark of Sino Biological Inc.



Profile of anti-FAS (CD95) reactivity on mouse splenocytes analyzed by flow cytometry. Cells should be Fc-blocked by treatment with Mouse BD Fc Block™ purified anti-CD16/CD32 mAb 2.4G2 (Cat. No. 553141) prior to staining, washed, then stained with APC Rabbit anti-FAS(CD95).

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to [www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html](http://www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html) for technical protocols.