

FAS / CD95 / APO-1 / TNFRSF6 Antibody (FITC), Rabbit MAb



Catalog Number: 50027-R004-F

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Mouse CD95 / APO-1 / TNFRSF6 / FAS protein (Catalog#50027-M08H)
Reagents:	FITC-conjugated Rabbit monoclonal antibody
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD95 / APO-1 / TNFRSF6 / FAS (rM CD95 / APO-1 / TNFRSF6 / FAS; Catalog#50027-M08H; NP_032013.2; Met1-Arg169) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.
Ig Type:	Rabbit IgG
Clone ID:	004
Specificity:	Mouse CD95 / APO-1 / TNFRSF6 / FAS
Concentration:	5 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.09% sodium azide
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Alternative Names:	FAS

APPLICATIONS

Applications:	FCM
----------------------	-----


RECOMMENDED CONCENTRATION

Please Note: Optimal concentrations/dilutions should be determined by the end user.

**FAS / CD95 / APO-1 / TNFRSF6
Antibody (FITC), Rabbit MAb**

Catalog Number: 50027-R004-F

EliteRmab® is a registered trademark of Sino Biological Inc.



Profile of anti-FAS (CD95) reactivity on mouse splenocytes analyzed by flow cytometry. Cells should be Fc-blocked by treatment with Mouse BD Fc Block™ purified anti-CD16/CD32 mAb 2.4G2 (Cat. No. 553141) prior to staining, washed, then stained with FITC Rabbit anti-FAS(CD95).

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.