

FAS / CD95 / APO-1 / TNFRSF6 Antibody, Rabbit MAb



Catalog Number: 50027-R006

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse CD95 / APO-1 / TNFRSF6 / FAS protein (Catalog#50027-M08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD95 / APO-1 / TNFRSF6 / FAS (rM CD95 / APO-1 / TNFRSF6 / FAS; Catalog#50027-M08H; NP_032013.2; Met 1-Arg 169).
Ig Type:	Rabbit IgG
Clone ID:	006
Specificity:	Mouse CD95 / APO-1 / TNFRSF6 / FAS
	No cross-reactivity in ELISA with Human FAS
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse FAS.

Please Note: Optimal concentrations/dilutions should be determined by the end user.