

FCGR2B / CD32 Antibody, Rabbit MAb



Catalog Number: 50030-R005

EliteRmab® is a registered trademark of Sino Biological Inc.

| GENERAL INFORMATION | |
|---------------------------|---|
| Immunogen: | 50030-M08H |
| Preparation | This antibody was obtained from a rabbit immunized with purified, recombinant Mouse FCGR2B / CD32 (rM FCGR2B; Catalog#50030-M08H; Met1-Arg217; NP_001070657.1). |
| Ig Type: | Rabbit IgG |
| Clone ID: | 005 |
| Specificity: | Mouse FCGR2B / CD32 |
| Formulation: | 0.2 µm filtered solution in PBS |
| Storage: | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles. |
| Alternative Names: | CD32A |
| APPLICATIONS | |
| Applications: | ELISA |
| RECOMMENDED CONCENTRATION | |
| ELISA | ELISA: 1:5000-1:10000 |

Please Note: Optimal concentrations/dilutions should be determined by the end user.