

CD32 / FCGR2B Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 50030-RP02

GENERAL INFORMATION

Immunogen:	50030-M08H
Preparation	Produced in rabbits immunized with purified, recombinant Mouse FCGR2B / CD32 extracellular domain (rM FCGR2B; Catalog#50030-M08H; Met1-Arg217; NP_001070657.1). FCGR2B specific IgG was purified by mouse FCGR2B affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse FCGR2B / CD32
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	CD32A

APPLICATIONS

Applications: ELISA

RECOMMENDED CONCENTRATION

ELISA ELISA: 1:1000-1:2000

Please Note: Optimal concentrations/dilutions should be determined by the end user.