

**CD62L / L-Selectin / SELL Antibody,  
Rabbit MAb**

Catalog Number: 50045-R023

EliteRmab® is a registered trademark of Sino Biological Inc.

**GENERAL INFORMATION**

<b>Immunogen:</b>	Recombinant Mouse CD62L / L-Selectin / SELL Protein (Catalog#50045-M08H)
<b>Preparation</b>	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CD62L / L-Selectin / SELL (rM CD62L / L-Selectin / SELL; Catalog#50045-M08H; NP_035476.1; Met1-Asn332).
<b>Ig Type:</b>	Rabbit IgG
<b>Clone ID:</b>	023
<b>Specificity:</b>	Mouse CD62L / L-Selectin / SELL
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

**APPLICATIONS**

<b>Applications:</b>	ELISA(Cap)
----------------------	------------

**RECOMMENDED CONCENTRATION**

<b>Sandwich ELISA (Capture Ab)</b>	ELISA(Cap): 1:250-1:2000 This antibody will detect Mouse CD62L / L-Selectin in ELISA pair set (Catalog: # SEK50045). In a sandwich ELISA, it can be used as capture antibody.
--	--

***Please Note: Optimal concentrations/dilutions should be determined by the end user.***